



NARKOMED 6000 TECHNICAL SERVICE BULLETIN # 6010

Product: Narkomed 6000

Date: October 1, 2001

Re: Phillips type back panel screws.

Reference Doc: N/A

Reason/Rational: Use of the traditional button head is not well suited for securing the Narkomed 6000 back panels. There are structural differences in the Narkomed 6000 that contribute to the difficulty of realignment and may lead to various mounting concerns. It is hard to apply sufficient torque to the button head to ensure a tight fit and could loosen, if started wrong may cross thread and sometimes break-off.

**Solution/
Procedure:** During routine PMC's inspect for loose or missing screws and replace as necessary with the new Phillips head screw, except for the bottom row of the lower back panel which shall continue to use the button head screws. During subsequent 2.06 upgrades or at any time the back panels are removed the button head screws shall be replaced with the Phillips head type screws, except for the bottom row of the lower back panel which shall continue to use the button head screws.

Devices affected: Phillips screws are now being used at the factory to secure all Narkomed 6000 back panels on shipments made after July 2001.

Timing/When: Immediately

Additional Info: Twenty Phillip head screws are required to secure the upper back panel. Thirteen Phillips screws are required to secure the lower back panel. The vapor park position panel and the small panel near the serial number are not affected by this change.

Costs: This screw is available at no charge if this TSB and the affected Narkomed 6000 serial numbers are referenced at the time the order is being placed.

**Ordering
Info/Parts:** P/N HW19003

Distribution: DrägerService, Authorized Service Organizations

**If you have any questions, please contact Technical Support by phone at
1-800-543-5047 or by fax at 1-215-721-5789**